

FOCAL DILATION WHERE IT MATTERS

Controlled
Growth


Consistent
Working
Range

Focal
Dilation

FOCAL DILATION WHERE IT MATTERS


Controlled Growth

Controlled Growth. Yangtze NC demonstrates a consistent controlled balloon growth across all diameters. Giving you confidence that you open each lesion and optimize each stent with the same efficiency.



Consistent Working Range


Consistent wide working range. Yangtze NC is designed with a wide working range enabling you to safely go to high pressure. It does so, regardless of diameters. Giving you confidence in safe balloon dilation.



Trademarks may be registered and are the property of their respective owners. Data on file at Stentys based on latest public available information.

Focal Dilation

Low compliance for focal dilation. By reducing the compliance, Yangtze NC is able to deliver optimal dilation forces to the area of most resistance to maximize lesion opening, while focusing all its power only where it matters for a safe opening of even very calcified lesions. Giving you confidence in safe lesion opening.



Trademarks may be registered and are the property of their respective owners. Data on file at Stentys based on latest public available information. May not be indicative of clinical performance.